



*“Thank you for your purchase and your trust in our product.”*

## - PRECAUTIONS AND WARNINGS –

Never use the facial derma roller on irritated or broken skin, or if you have active skin conditions such as severe acne.

If you have specific skin concerns, consult a healthcare professional or dermatologist.

For hygiene reasons, do not share your facial derma roller.

It is generally recommended to use the facial derma roller 1–2 times per week, depending on your skin’s tolerance.

**KEEP OUT OF REACH OF CHILDREN**

# - HOW TO USE THE FACIAL DERMA ROLLER -

## 1. PREPARATION :

Disinfect your facial derma roller by soaking it for a few minutes in 90% medical isopropyl alcohol.

Rinse it and allow it to air-dry completely.

In the meantime, cleanse your face thoroughly.

## 2. APPLICATION :

Use the device preferably in the evening, as your skin may appear slightly red afterwards.

Apply light pressure and gently roll the derma roller vertically, horizontally, and diagonally, avoiding sensitive areas such as the eyes and lips.



Make 2–4 passes in each direction and on each area.

Do not exceed 2 minutes of treatment time.

### 3. AFTERCARE :

Once you're done, you can apply a hydrating or soothing serum to support skin comfort. Check the ingredient list to avoid any potential irritation.

Avoid using harsh products such as exfoliants or retinoids immediately after treatment.

Stay out of direct sunlight and use sunscreen during the days following the session.

### 4. MAINTENANCE :

After each use, disinfect your facial derma roller again and let it air-dry.

Store it in its case with the lid closed.

For optimal hygiene and performance, replace the roller heads every 3 to 6 months.

*“Feel free to leave us a review or share your feedback. ”*